

## BIOMAG-I

# Clinical outcomes of the third-generation resorbable magnesium scaffold for coronary artery lesions: three-year results of the BIOMAG-I<sup>1</sup>

## Conclusions

- The low TLF rate at three years – particularly with only one event occurring beyond the scaffold resorption period – is highly encouraging.<sup>1</sup>
- Freedom from cardiac death, target vessel myocardial infarction, and any device thrombosis up to three years attests to **excellent device performance**.<sup>1</sup>
- The favorable 3-year outcomes support renewed interest in bioresorbable scaffolds as a viable therapeutic option that combines **temporary mechanical support with excellent long-term safety and efficacy**.<sup>1</sup>

## Study objective

Assessment of angiographic, clinical and safety performance of the **sirolimus eluting magnesium scaffold Freesolve** in patients with de novo coronary artery lesions

## Primary endpoint

- In-Scaffold LLL at 6 months

## Secondary endpoints at 12 months

- Angiographic in-scaffold and in-segment LLL
- Binary restenosis
- Diameter stenosis
- IVUS/OCT descriptive analysis

## Clinical endpoints

- Target Lesion Failure\*
- Definite or Probable Scaffold Thrombosis


## Patient characteristics

	n = 116	%
Age, years	61.0 ± 9	
Male	90	77.6%
Hypertension	86	74.1%
Hypercholesterolemia	72	62.1%
Diabetes	32	27.6%
History of smoking	75	64.7%
History of myocardial infarction	39	33.6%
<b>NSTEMI</b>	<b>24</b>	<b>20.7%</b>

## Lesion characteristics

	n
Lesion length (mm)	12.3 ± 5.1
Reference vessel diameter (mm)	2.72 ± 0.46
<b>AHA/ACC lesion class B2/C</b>	<b>90</b>
Side branch involvement	25


Lesion characteristics are estimated by core lab.




## Lesion location

	n	
LAD	53	45.3%
LCx	22	18.8%
RCA	40	34.2%
Ramus intermedius	2	1.7%

## In-Scaffold Late Lumen Loss at 12 months<sup>2</sup> (n = 100)



## Target Lesion Failure



## Serial OCT Analysis at 12 months



## Coordinating investigator

Michael Haude, MD, PhD, Rheinland Klinikum, Neuss, Germany

\*TLF is defined as Composite of Cardiac Death, TV-MI, CD-TLR [Kaplan-Meier estimate]; \*\*peri-procedural target vessel MI according to SCAI definition and non-peri-procedural target vessel MI according to Universal MI Definition; <sup>†</sup>driven by three clinically-driven target lesion revascularization; <sup>‡</sup>Definition of malapposition: if the distance between outer contour of the strut and vessel wall is more than the individual strut thickness; <sup>§</sup>p < 0.05 for 12-month vs post-procedure. 1. Haude et Al. EuroIntervention 2025;21:e1-e3. 2. Haude et Al. EuroIntervention 2023;19:e414-e422. 3. Haude et Al. Eur Heart J 2016;37:2701-2709.

All endpoint related events have been adjudicated by an independent clinical event committee. BIOMAG-I and BIOSOLVE-II are based on Kaplan-Meier failure estimate analysis including censored observations.

Freesolve (DREAMS 3G) and Magmaris (DREAMS 2G) are trademarks or registered trademarks of the BIOTRONIK Group of Companies. Commercially available in selected countries.