



Head-to-head comparison of **Passeo[®]-18 Lux[®]** DCB to IN.PACT Admiral DCB at 12 Months¹

Conclusions

- The BIOPACT RCT demonstrated the noninferiority of the low-profile² **Passeo-18 Lux** DCB compared to IN.PACT Admiral, with excellent safety and efficacy results at 12 months:
 - Freedom from CD-TLR³: **Passeo-18 Lux**, 97.2% vs IN.PACT Admiral, 97.0% (p = 0.0002)
 - Primary safety endpoint⁴: **Passeo-18 Lux**, 95.7% vs IN.PACT Admiral, 96.3% (p = 0.0008)
 - Primary patency⁵: **Passeo-18 Lux**, 88.2% vs IN.PACT Admiral, 89.1% (p = 1.00)

Study design

Physician-initiated, prospective, multicenter, single-blinded, 1:1 randomized controlled non-inferiority trial comparing the safety and efficacy of the Passeo-18 Lux drug-coated balloon versus the IN.PACT Admiral drug-coated balloon for treatment of stenotic, restenotic or occlusive lesions of the femoropopliteal arteries.

Endpoints

Primary Endpoints


- Efficacy: 12-month freedom from CD-TLR²
- Safety: composite of freedom from device- and procedure-related death through 30 days, freedom from major target limb amputation and clinically driven target vessel revascularization through 12 months

Secondary Endpoints

- 6- and 12-month freedom from MAE (all-cause death, clinically driven target vessel revascularization, major target limb amputation, or thrombosis at target lesion)
- 6- and 12-month primary patency
- 6- and 12-month freedom from all-cause mortality
- 6- and 12-month freedom from major target limb amputation

Patient characteristics

Patient characteristics	Passeo-18 Lux n = 152	IN.PACT Admiral n = 149
Age ± SD, years (range ± SD)	69 ± 8 [47-87]	67 ± 9 [44-90]
Diabetes mellitus	39 [25.7%]	44 [29.5%]
Hypertension	109 [71.7%]	118 [79.2%]
Hypercholesterolemia	110 [72.4%]	116 [77.9%]
Smoking	58 [38.2%]	74 [49.7%]
Previous arterial intervention	73 [48%]	55 [36.9%]
Previous coronary intervention	42 [27.6%]	37 [24.8%]
Rutherford Classification		
Class 2 Moderate	39 [25.8%]	35 [23.5%]
Class 3 Severe	101 [66.9%]	103 [69.1%]
Class 4 Ischemic Rest Pain	11 [7.3%]	11 [7.4%]



Lesion characteristics


Lesion characteristics	Passeo-18 Lux	IN.PACT Admiral
Reference vessel diameter [mm]*	5.3 ± 0.7 [4-7]	5.4 ± 0.7 [4-7]
Lesion length (mm)*	74 ± 49.4 [4-180]	65.6 ± 39 [5-180]
% Stenosis (mm)*	84.1 ± 9.7 [50-99]	84.8 ± 9.4 [70-99]
Total occlusion	26 [17.1%]	31 [20.8%]
Calcification [moderate-severe]	121 [79.6%]	106 [71.1%]
Lesion location		
Proximal SFA	16 [10.5%]	15 [10.1%]
Mid SFA	71 [46.7%]	76 [51.0%]
Distal SFA	68 [44.7%]	54 [36.2%]
P1	16 [10.5%]	21 [14.1%]

Procedure details

Procedure details	Passeo-18 Lux	IN.PACT Admiral
Pre-dilatation	152 [100%]	149 [100%]
Bailout stenting	16 [10.5%]	17 [11.4%]
Postdilatation	23 [15.1%]	20 [13.4%]


* Data shown as mean ± SD (range)

Freedom from CD-TLR²




Number at risk	0	90	180	270	335	365	395
IN.PACT Admiral	149	141	139	133	132	104	76
Passeo-18 Lux	152	149	145	141	137	114	78

Primary patency⁵ (core lab adjudicated)



Number at risk	0	90	180	270	335	365	395
IN.PACT Admiral	149	135	125	109	108	101	99
Passeo-18 Lux	152	137	131	113	112	109	105


Primary safety endpoint⁴



Number at risk	0	90	180	270	335	365	395
IN.PACT Admiral	149	141	139	132	131	104	76
Passeo-18 Lux	152	148	144	140	135	113	78




Event	Passeo-18 Lux	IN.PACT Admiral
Death related to procedure/device	0	0
Target vessel revascularization	6	5
Major amputation	0	0

Survival



Number at risk	0	90	180	270	335	365	395
IN.PACT Admiral	149	141	139	133	132	104	76
Passeo-18 Lux	152	149	145	141	137	114	78

Passeo-18 Lux key outcomes



Principal investigator

Koen Deloose, Dendermonde, Belgium

CD-TLR: Clinically Driven Target Lesion Revascularization, DCB: Drug-Coated Balloon, MAE: Major Adverse Events, P1: Popliteal Segment 1, SFA: Superficial Femoral Artery.

1. BIOPACT RCT: a third head to head DCB RCT. Deloose K. Presented at Paris Vascular Insights 2022; 2. BIOTRONIK data on file; 3. Freedom CD-TLR: Freedom from any reintervention at the target lesion due to symptoms, drop of ankle-brachial index (ABI) > 20% or ABI > 0.15 compared to postprocedural ABI; 4. Primary safety endpoint: composite of freedom from device- and procedure-related death through 30 days, freedom from major target limb amputation and clinically driven target vessel revascularization through 12 months; 5. Primary patency: composite of freedom from CD-TLR and binary restenosis (restenosis defined as duplex ultrasound [DUS] peak systolic velocity ratio ≤ 2.4 or $\leq 50\%$ stenosis as assessed by an independent DUS core lab at 12-month follow-up).

The **Passeo-18® Lux® DCB** with its **Lux®** coating is part of the **Lux®** family of Paclitaxel-coated balloons from BIOTRONIK.

Passeo and Lux are trademarks or registered trademarks of the BIOTRONIK Group of Companies. IN.PACT Admiral is a trademark or registered trademark of Medtronic.