





27-JAN-25


















Follow-Up Protocol		
2023 HRS/APHRS/LAHRS guideline on cardiac physiologic pacing for the avoidance of heart failure	and mitigation	
In patients who have had a CSP device implanted, a 12-lead ECG, including long strips during threshold testing, can help to ensure and optimize maximal conduction system capture. The tracing should be evaluated to determine capture thresholds, LBBB correction when pertinent, and type of capture (selective vs nonselective conduction system capture). The paced QRS duration and morphology should be compared to prior readings and used as a comparison point for future follow-up ¹		
1. Chung MK, Patton KK, Lau CP, et al. 2023 HRS/APHRS/LAHRS guideline on cardiac physiologic pacing for the avoidance and mitigation of heart failure. <i>Heart Rhythm</i> . 2023;20(9):e17-e91. doi:10.1016/i.hrthm.2023.03.1538		
17 JANUARY, 2025	BIOTRONIK excellence for life	








































Not BIOTRONIK Specific Table 1: Programming Recommendations for LBBAP				
Parameter	Recommendation HBP	Recommendation LBBAP		
Pacing mode	Single-chamber device: VVI Dual-chamber or CRT device: • HBP lead in a ventricular port: DDD(R), DDI or managed ventricular pacing mode • HBP lead in atrial port (chronic AF) with backup ventricular lead: DDD(R), DDI(R), or DVI(R) if available	Single-chamber device: VVI Dual-chamber or CRT device: • LBBAP lead in a ventricular port: DDD(R), DDI or managed ventricular pacing mode • LBBAP lead in atrial port (chronic AF) with backup ventricular or ICD lead: DDI(R), or DVI(R) if available		
Pacing polarity	Unipolar (better visibility or the pacing spike to avoid confounding with intrinsic rhythm, lower capture thresholds) Bipolar (lower current drain due to higher impedance)	Bipolar (lower current drain due to higher impedance; anodal capture may narrow the QRS) Unipolar if anodal capture is not desirable		
Sensing vector	Bipolar (unipolar can be tried if low sensing amplitude or P wave/HB potential oversensing)	Bipolar		
Sensitivity	HBP lead connected to atrial channel: set to the maximum value (minimum sensitivity), as ventricular sensing is provided by the backup ventricular lead HBP lead connected to RV channel: adjust the level to ensure ventricular sensing, yet avoid oversensing of atrial or HB potentials.	Usually, not an issue as R waves are of high amplitude		
Output voltage	2 x threshold voltage Fixed safety margin, e.g. 1V above the threshold, in non-dependent patients	2 x threshold voltage		
Impulse duration	$0.4\ {\rm ms}$ (1.0 ms if high capture threshold). 0.2-0.4 ms may be programmed according to chronaxle	0.4 ms (capture threshold is rarely an issue)		
Automatic capture control algorithms	Deactivate, monitoring only (may be inaccurate or impossible to measure, especially if the HBP lead is connected to the atrial port), or activate only once the accuracy has been confirmed in the patient	Set to monitor or automatic once the accuracy has been confirmed in the patient		
AV delay	HBP lead in ventricular port: Subtract HV interval (e.g. 40 ms) from desired AV interval	LBBAP lead in ventricular port: Subtract LBB-V interval (e.g. 2) ms) from desired AV interval		

Programming LBBAP Per Physician Preference (continued) Not BIOTRONIK Specific

AV delay	HBP lead in ventricular port: Subtract HV interval (e.g. 40 ms) from desired AV interval HBP lead in atrial port with backup ventricular pacing: AV delay >His pace- RVSS interval (e.g. 150 ms) HBP lead in atrial port with HOT-CRT: optimize AV interval based on QRS	LBBAP lead in ventricular port: Subtract LBB-V interval (e.g. 20 ms) from desired AV interval LBBAP lead in atrial port with backup ventricular pacing: AV delay >LBBAP- RVS interval (e.g. 150 ms) LBBAP lead in atrial port with LOT-CRT: optimize AV interval based on QRS
	narrowing, or program empirically to 60% of the His pace – RV sense interval (usually 40-60 ms) $^{\rm 8}$	narrowing
VV delay (CSP lead connected to LV port)	With backup RV pacing: program maximum LV channel pre-excitation (e.g. 80 ms) In case fusion with RV pacing is desirable (e.g. in case of uncorrected RBBB): program LV channel pre-excitation 30-60 ms, optimized by surface ECG	With backup RV pacing (e.g. with ICD lead): program maximum LV channel pre-excitation (e.g. 80 ms)
Ventricular safety pacing	Deactivate if the HBP lead is connected to the atrial port with an RV back-up lead, after having verified absence of crosstalk	Deactivate if the LBBAP lead is connected to the atrial port with an RV back- up lead (e.g. with an ICD or in case of LOT-CRT), after having verified absence of crosstalk
Automatic sensing control algorithms	Deactivate (IP wave oversensing and HB sensing (may lead to asystole!)	Can be left on
Sensing if CSP lead connected to LV port	Deactivate (BIOTRONIK, Boston-Scientific)	Deactivate (BIOTRONIK, Boston-Scientific)
AV and VV optimization algorithms	Deactivate	Deactivate
Ventricular triggered pacing (ventricular sense response, etc.)	Deactivate	Deactivate
	system pacing; HB = His bundle pacing; HV = His-ventricle; HOT-CRT = His-optimized CR eft bundle branch pacing optimized CRT; LV = left ventricular; RBBB = right bundle branch ricular.	
30 JANUARY, 2025		BIOTRONIK excellence for life

Quiz Questions
1. Which of the following ECG characteristics best describes ns-LBBP?
a. Fusion between S-LBBP and LVSP, clear R' in V1 and narrow QRS
b. Clear R' in V6
c. Latency interval between pacing spike and QRS onset (not always clear)
d. None of the above
2. Which of the following ECG characteristics best describes s-LBBP?
a. Fusion between S-LBBP and LVSP
b. Large R' in V1, narrow QRS and Latency interval (separation) between pacing spike and QRS onset (not always clear)
c. Wide QRS with notching
d. None of the above
31 JANUARY, 2025 SIGTRONIK

