



Amvia Edge HF-T QP  
Home Monitoring

1234567890

**BIOTRONIK**  
Made in Germany

0AE-DDDRV

RA IS-1  
LV IS4 LLLL  
RV IS-1

## Rate Adaptive Sensors: Past & Present

David Hayes, M.D. – Chief Medical Officer  
Doug Finch – Sr. Advanced Product Program Manager




**BIOTRONIK**  
excellence for life

1

## Characteristics of Ideal Sensor

- Appropriate speed of acceleration/deceleration
- Proportionality to workload
- Sensitive to exercise and non-exercise requirements, i.e. postural changes, emotional states, circadian variations
- Specificity, i.e. avoid false sensor activation by non-physiological needs


2 | FEBRUARY 20, 2024 |



**BIOTRONIK**  
excellence for life

2

## Open Loop vs. Closed Loop Sensors




3 | FEBRUARY 20, 2024 |



3

## Schematic of Physiologic Drivers



4 | FEBRUARY 20, 2024 |



4

## Sensors

### Historical sensors

- *Piezoelectric crystal* — — — — — ▶
- Temperature
- QT interval
- O2 Saturation
- Pre-ejection interval
- Ventricular depolarization gradient
- *Peak endocardial acceleration*

### Currently available

- Accelerometer
- Minute ventilation
- Myocardial Contractility (CLS)



















5


## Properties of a Normal Sinus Node

1. **Increases HR with exercise**
  - Proportionate and physiologically controlled
2. **Acute mental stress**
3. **P-P variability**
4. **Diurnal variability**
  - “normal” HR distribution/24 Hrs.
5. **Circadian variability**
  - “normal” HR distribution during sleep
6. **Compensatory response**
  - Changes in MABP due to vasodilation

6

### Comparison of Sensor Ability to Mimic a Normal Sinus Node

	Accelerometer	MV	CLS
Increases HR with exercise			
Acute mental stress			
P-P variability			
Diurnal variability			
Circadian variability			
Compensatory response			

12-Feb-25 

7

## What is CLS & How Does It Work?


- **Primary Sensor**
  - Takes measurements directly from the myocardial tissue
- **Built on a beat-to-beat architecture<sup>1</sup>**
  - Rapidly evaluates every beat and is able to respond on the next beat
- **Automated with machine learning<sup>2</sup>**
  - Automatically learns the patient
  - Automatically recalibrates to systemic variations, e.g. changes in hydration, medication, etc.

*Takes measurements of myocardial contractility (not motion or breathing). Why is this important? HR & BP should be correlated and are both managed by the autonomic nervous system*

*No sudden change criteria, resulting in immediate response to change. Why is this important? That's how the autonomic nervous system manages HR*


*No need for optimization, fully automated. Why is this important? CLS incorporates the device into the autonomic nervous system.*

<sup>1</sup>, Malinowski et al. PACE 1998;21:11  
<sup>2</sup>. Zechhi P et al. Prog Biomed Res 2000;5:2

12-Feb-25 

8

## CLS – Impedance Principle




9

9

## Impedance Changes into Rate

Increase in area differential

Physiologic increase in paced rate





10

10

# Closed Loop Stimulation (CLS)

- CLS patients showed a 75% reduction in orthostatic hypotension events <sup>1</sup>
- Decreases in blood pressure can be modulated by CLS <sup>1, 2, 3</sup>
- CLS achieves a reduction in reoccurring syncope <sup>7, 8</sup>
- CLS precisely correlates heart rate and blood pressure <sup>11</sup>
- CLS responds to acute mental stress <sup>6</sup>
- CLS significantly reduces atrial burden <sup>4, 5</sup>
- CLS reduced frequency of AT episodes less than one hour <sup>4</sup>
- CLS achieves heart variability that is analogous to a healthy sinus node <sup>9, 10</sup>



1. Abi-Samra FM et al. *Europace*. 2013, 15(6).
2. Quaglione R et al. *Europace*. 2005, 7(4).
3. Palmisano P et al. *Europace*. 2017, Epub.
4. Puglisi A et al. *PACE*. 2008, 31(11).
5. Puglisi A et al. *Eur Heart J*. 2003, 24(21).
6. Chandiramani S et al. *PACE*. 2007, 30(8).
7. Baroni-Esquivias, G. *ACC* 2017, 70(14).
8. Occhetta E et al. *Europace*. 2004, 6(6).
9. Quaglione R et al. *PACE*. 2009, 33(3).
10. Malinowski K. *PACE*. 1998, 21(11).
11. Zecchi P et al. *Prog Biomed Res*. 2000 5(2).
12. Tse HF et al. *JACC*. 2005, 46(12).
13. Wojciechowski D et al., *Prog Biomed Res*. 2001, 6.
14. Novak M et al. *Prog in Biomed Res*. 1998, June.

11

# Closed Loop Stimulation

## VICRA Study

- CLS is the only pacemaker algorithm to respond to physiologic demands and acute mental stress on a beat-to-beat basis
- CLS precisely correlates heart rate and blood pressure

- 1 CLS delivers gradual rate decrease based on reduction of myocardial contractility.<sup>1</sup>
- 2 Only CLS provides physiological rate adaptation in response to acute mental stress.<sup>1</sup>
- 3 CLS provides rate adaptation during hemodynamic changes induced by Isoproterenol infusion.<sup>1</sup>




12 1. Zecchi P et al. *Prog Biomed Res*. 2000 5(2).




12

## CLS vs. VVIR vs. VVI

Mean HR Increase



Mean Pacing (%)



\*Tests of within-subject contrasts

Proietti et al: PACE 2012;00:1-9

13 | FEBRUARY 20, 2024



13

## Case Study 1



You are performing a follow up on a 75-year-old patient implanted with an EDORA 8 DR-T pacemaker implanted 3 years ago. She is feeling well and is not expected to require any programming changes today. When asked, she notes that she is tired and really never has much energy but feels that is pretty normal for her.


She is Atrial paced 99% and Ventricular paced 100%. There is nothing remarkable from her previous follow-ups.

14 | FEBRUARY 20, 2024



14

## Case Study 1



Programmed:

- DDD 60 bpm /130 bpm
- AV Delay 180 ms-140ms (nominal)

Which of the following would be next most helpful in assessing her fatigue?

1. Ventricular pacing safety margin
2. A & V histograms
3. Atrial arrhythmia burden
4. Ventricular lead impedanc

## Case Study 1



Diagnostics:

- Atrial histogram - Blunted – 99% at Base Rate
- Ventricular histogram - Blunted – 100% V Paced



## Case Study 1




### Sensing Test:

P & R waves within normal range

Which of the following parameters necessitates alternate programming?

1. Atrial sensitivity
2. Post-atrial ventricular blanking period
3. AV interval
4. Ventricular pulse-width

## Case Study 1



### Final Programming

Two crucial programming changes were made for this patient considering her blunted rate histograms; 99% A & 100% V pacing; and intact intrinsic P & R waves.


1. Mode changed to DDDR (DDD-CLS)
  - Either CLS or accelerometer could be used but CLS was chosen given physiologic advantages of CLS
  - Consider the patient's age (220-age X 80%) when deciding Max sensor rate
2. AV Delay optimization: either I-OPT or AV Hysteresis

## Case Study 2



You are called to perform a routine follow-up on a 65-year-old patient with a dual chamber ICD. The patient's device has not been interrogated for over a year and the patient complains of increasing fatigue during normal daily activities despite optimal medical therapy and no apparent decompensation. No programming changes are expected to be made today.

## Case Study 2



What do we learn from the sensing test?

## Case Study 2



With a blunted histogram, what programming options should be considered?

21 | FEBRUARY 20, 2024 |



21

## Case Study 2



### Final programming options

- DDD-CLS to improve fatigue
- I-Opt to maintain low RV pacing percentage

22 | FEBRUARY 20, 2024 |



22

## Case Study 2



- Turning CLS on by changing the mode to DDD-CLS
- Turn on I-Opt through the AV delay screen and choosing I-Opt from the AV hysteresis mode


23 | FEBRUARY 20, 2024 |



23

## Simply Program CLS on by Selecting DDD-CLS Mode

Programming: One button to program CLS "On"




24 | FEBRUARY 20, 2024 |



24

## Set the Basic Rate and the Max CLS Rate, and CLS Will Automatically Optimize to the Patient

Programming: Easier to program than an accelerometer





25 | FEBRUARY 20, 2024 |



25

## How Do I Know This Patient is Benefiting From CLS?

- Atrial rate histogram demonstrates patient was achieving pacing rates from 60 to 130, driven by CLS
- Sensor histogram shows the accelerometer would have paced at the lower rate limit of 60 bpm 100% of the time because no motion was detected
- A quick look at the activity report would verify that the patient is indeed sedentary




26 | FEBRUARY 20, 2024 |




26

# Understanding Rate Distribution

Rate Distributions for CLS Dynamics Settings



Atrial Rate Histograms at different CLS Response settings



27

# Expert Options

The screenshot shows the ICD programming software interface. At the top, there are ECG waveforms for A, V, FF, A, and V channels. Below the waveforms, a status bar indicates "Programming was successful." The main menu includes Tachycardia, Bradycardia, Home Monitoring, Diagnostics, Patient, and MRI. The "CLS" (Cardiac Loop Sensing) expert options dialog box is open, showing the following settings:
 

- Mode: DDD-CLS
- Basic rate [bpm]: 50
- Pulse amplitude [V]: 3.5
- Max. CLS rate [bpm]: 120
- CLS response: Medium
- CLS resting rate control [bpm]: +20
- Vp required: No

 The interface also shows various system status indicators on the left (e.g., 4G, LiveControl, USB) and a sidebar on the right with options like ICD therapy, RF telemetry, Follow-up, Parameters, Tests, Recordings, Diagnostics, Status, LiveSupport, More, Preferences, and End.

28

## CLS Response Programming

### Expert options

The screenshot shows the Biotronik programmer interface for an Ilivia 7 DRT device. The main window is in the 'Home Monitoring' tab, showing 'DDD-CLS' mode with a basic rate of 50 bpm. A 'CLS response' dialog box is open, allowing selection of a response level: Very low, Low, Medium (highlighted), High, or Very high. The interface also shows various monitoring and therapy settings on the right side.

29 | FEBRUARY 20, 2024 |

BIOTRONIK  
excellence for life

29

## Suggested Programming Modifications for Syncope Patients

### Resting rate control

- Resting rate control prevents inappropriate rate excursions during orthostatic changes or response to non-physical stimuli
- Resting rate control should be turned OFF only for dysautonomia patients:
  - Neurocardiogenic / cardio-inhibitory syncope
  - Vasovagal syncope
  - POTS
  - Autoimmune dysautonomia

The screenshot shows the Biotronik programmer interface with the 'CLS resting rate control' dialog box open. The dialog box allows selection of a resting rate control level: OFF, +10, +20, +30, +40, or +50 bpm. The background shows the 'Home Monitoring' tab with 'DDD-CLS' mode and a basic rate of 50 bpm.

30 | FEBRUARY 20, 2024 |

BIOTRONIK  
excellence for life

30

## Rate-Adaptive Sensors: Past & Present

- Multiple sensors have been tried but only 3 are currently market available
- Sensors can be 'open-loop' or 'closed-loop' and clinically it is important to understand the difference
- CLS is the only 'closed loop' sensor available and has been shown to be advantageous for multiple clinical indications
- While programming CLS can be as simple as 'on' there are some more advanced programming options that can help optimize rate adaptation for an individual patient