

# Combining Home Monitoring Temporal Trends from Implanted Defibrillators and Baseline Patient Risk Profile to Predict Heart Failure Hospitalizations (SELENE HF)

D'ONOFRIO A ET AL.,  
EUROPACE 2021

## Study Design

- International, multi-center, prospective, observational, event-driven cohort study
- 34 centers in Italy and Spain
- 918 ICD (dual-chamber, DX and CRT-D) patients with NYHA Class II-III, LVEF  $\leq$  35%
- To develop and validate an algorithm to predict heart failure (HF) hospitalization based on 7 Home Monitoring parameters and a baseline risk stratifier (based on the Seattle Heart Failure Model)
- Home Monitoring parameters: Mean heart rate, mean heart rate at rest, premature ventricular contractions (PVC), atrial burden, heart rate variability (HRV), patient activity, and thoracic impedance (TI)
- Primary endpoint: First HF-related hospitalization
- Collection of  $\geq$  50 primary endpoint events
- Post-hoc randomization into 2 cohorts for algorithm development and validation

## Main Results

The algorithm predicted HF hospitalizations early, with high sensitivity and with a low false alert rate.

Endpoint	Sensitivity (%)	Alerting time (days)	False alert rate (ppy)
First post-implant HF hospitalization	65.5 (45.7-82.1)	42 (21-89)	0.69 (0.64-0.74)

Table 1: Main performance parameters of the predicting algorithm  
(Numbers in brackets: confidence intervals, except for alerting time (IQR))

### 66% sensitivity

2 out of 3 hospitalizations were predicted

### 42 days alerting time

may allow proactive care to possibly prevent hospitalizations

### Only 0.7 false alerts per patient year


may keep workload low

## Clinical Relevance

- HF is associated with poor prognosis and high hospitalization rates. Recurrent hospital admission due to HF result in a gradual worsening of the health status of patients and constitutes a considerable healthcare burden.
- Early prevention of HF decompensation is a key strategy to improve patient outcomes. A CIED-based algorithm that predicts impending HF hospitalization could help reduce hospitalization among high-risk HF patients.

# The Newly Designed Algorithm is Early and Highly Effective in Predicting HF Hospitalizations

## Patient case: Prediction of a HF Hospitalization Event



Trends of Home Monitoring variables and the predicting index in an 82-year-old man preceding a hospitalization event.

Relevant changes in trends of **multiple Home Monitoring parameters** contribute to a **single predicting index**.

The index increases over time, ultimately leading to an alert.

**The alerting time of 42 days** prior to hospitalization would have allowed proactive care and possibly would have prevented the hospitalization.

### Main drivers of the index increase:

- Mean heart rate (24h HR)
- Mean heart rate at rest (Night HR)
- Patient activity (24h Activity)
- Premature ventricular contractions (PVC/day)
- Thoracic impedance (TI)

### The predicting index appears in the box:

- Day of alert
- Day of HF hospitalization
- Threshold of predicting index
- In-alert state (period of increased risk of HF hospitalization)

## What the Authors Say:

“ Our algorithm generated alerts at a median of 42 to 61 days before events, allowing sufficient time for patient contact, investigations, and preventive measures to reduce hospitalizations. ”

Source: D'Onofrio et al. Combining home monitoring temporal trends from implanted defibrillators and baseline patient risk profile to predict heart failure hospitalizations: results from the SELENE HF study. *Europace* 2021, doi: 10.1093/europace/euab170.