Pantera LEO

English	Fast-Exchange PTCA catheter
Deutsch	Fast-Exchange-PTCA-Katheter
Français	Cathéter pour ACTP à échange rapide
Italiano	Catetere PTCA a scambio rapido
Español	Catéter de intercambio rápido para ACTP
Български	PTCA катетър за бърза смяна
Hrvatski	PTCA kateter za brzu izmjenu
Česky	Katetr PTCA pro rychlou výměnu
Dansk	PTCA-kateter med fast exchange
Nederlands	Snel verwisselbare PTCA-katheter
Eesti	Kiirelt vahetatav PTCA-kateeter
Suomi	Nopeasti vaihdettava PTCA-katetri
Ελληνικά	Καθετήρας PTCA ταχείας εναλλαγής
Magyar	Gyorsan cserélhető (Fast-change) PTCA katéter
Latviešu	Ātri nomaināms PTCA katetrs
Lietuvių	Greitai keičiamas PTCA kateteris
Norsk	PTCA-kateter for rask utskiftning
Polski	Szybkowymienny cewnik do PTCA
Português	Cateter de troca rápida para PTCA
Română	Cateter PTCA cu schimbare rapidă
Русский	Быстрозаменяемый ЧТКА-катетер
Slovenčina	Katéter na PTCA s funkcioi rýchlej výmeny
Slovenščina	Kateter PTCA za hitro zamenjavo
Svenska	PTCA-kateter för snabbt byte
Türkçe	Hızlı Değişirme Özelliğine Sahip PTCA kateteri

Instructions for Use	7
Gebrauchsanweisung	8
Mode d'emploi	9
Manuale tecnico di istruzioni	11
Instrucciones de uso	12
Ръководство за употреба	14
Priručnik za uporabu	15
Návod k použití	17
Brugsanvisning	18
Technische handleiding	20
Kasutusjuhend	21
Käyttöohje	22
Οδηγίες χρήσης	24
Használati utasítás	25
Lietošanas norādījumi	27
Naudojimo instrukcija	28
Bruksanvisning	30
Instrukcja obsługi	31
Manual técnico	32
Instrucțiuni de utilizare	34
Инструкция по применению	35
Návod na použitie	37
Navodila za uporabo	38
Bruksanvisning	40
Kullanım Talimatları	41

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

BIOTRONIK
excellence for life

Compliance Chart/ Elastizitätstabelle/ Tableau de compliance/ Tabella di compliance/ Tabla de distensibilidad/ Таблица за разтегливост/ Tablica skladnosti/ Tabulka poddajnosti/ Kompianskort/ Compliantietabel/Vastavustabel/ Compliance Chart (venyyvystaulukko)/ Διόγραφο ενδοτικότητας/ Kompiansz táblázat/ Tabula atilstības noteikšana/ Atitikties lentelė/ Kompatibilitetstabel/ Tabla zgodności/ Tabela de complacência/ Tabel de conformitate/ Таблица растяжимости/ Tabulka poddajnosti/ Preglednica skladnosti/ Elasticitetstabel/ Esneklik Tablosu

Inflation Pressure (IP)			Balloon Diameter (mm) (BD)										
	atm	(kPa)	2.0	2.25	2.5	2.75	3.0	3.25	3.5	3.75	4.0	4.5	5.0
	6	(608)	1.93	2.16	2.41	2.65	2.89	3.12	3.36	3.57	3.85	4.38	4.86
	8	(811)	1.95	2.18	2.43	2.68	2.92	3.15	3.40	3.61	3.89	4.41	4.90
	10	(1013)	1.96	2.21	2.45	2.70	2.94	3.19	3.43	3.66	3.92	4.44	4.93
	12	(1216)	1.98	2.23	2.48	2.73	2.97	3.22	3.47	3.70	3.96	4.47	4.97
NP	14	(1419)	2.00	2.25	2.50	2.75	3.00	3.25	3.50	3.75	4.00	4.50	5.00
	15	(1520)	2.01	2.26	2.51	2.76	3.01	3.27	3.52	3.77	4.02	4.52	5.02
	16	(1621)	2.02	2.27	2.52	2.77	3.03	3.28	3.54	3.80	4.04	4.53	5.04
	17	(1723)	2.03	2.28	2.54	2.79	3.04	3.30	3.55	3.82	4.06	4.55	5.05
RBP	18	(1824)	2.04	2.30	2.55	2.80	3.06	3.32	3.57	3.84	4.08	4.56	5.07
	19	(1925)	2.05	2.31	2.56	2.81	3.07	3.33	3.59	3.87	4.10	4.58	5.09
RBP	20	(2027)	2.05	2.32	2.57	2.83	3.09	3.35	3.61	3.89	4.12	4.60	5.11
	21	(2128)	2.06	2.33	2.58	2.84	3.10	3.37	3.62	3.92	4.14	4.61	5.12
	22	(2229)	2.07	2.34	2.60	2.85	3.11	3.38	3.64	3.94	4.16	4.63	5.14
NP	In vitro testing has shown that the balloons will reach their nominal size at given Nominal Pressure.												
RBP	In vitro testing has shown that with 95% confidence, 99.9% of the balloons will not burst at or below Rated Burst Pressure. Do not exceed RBP.												
NOTE	Diameters at inflation pressures above RBP are provided for information only. At the pressures above RBP shown in this table, 99% of the balloons with 95% confidence will not burst.												

Sizes/ Größen/ Tailles/ Dimensioni/ Tamaños/ Размеры/ Veličine/ Velikosti/ Størrelser/ Afmetingen/ Suurused/ Koot/ Μεγέθη/ Méretek/ Izmēri/ Dydziai/ Størrelser/ Rozmiary/ Tamanhos/ Dimensiuni/ Размеры/ Vetkosti/ Velikosti/ Storlekar/ Boyutlar

Balloon Length (mm) (BL)	Balloon Diameter (mm) (BD)											
	2.0	2.25	2.5	2.75	3.0	3.25	3.5	3.75	4.0	4.5	5.0	
8	X	X	X	X	X	X	X	X	X	X	X	X
12	X	X	X	X	X	X	X	X	X	X	X	X
15	X	X	X	X	X	X	X	X	X	X	X	X
20	X	X	X	X	X	X	X	X	X	X	X	X
30	X	X	X	X	X	X	X	X	X	X	X	X

Abbreviations/ Abkürzungen/ Abréviations/ Abbreviazioni/ Abreviaturas/ Абревиатури/ Kratica/ Zkratky/ Forkortelser/ Afkorting/ Lühendid/ Lyhenteet/ Συμτομογραφίες/ Rövidítések/ Saisinājumi/ Santrumpos/ Forkortelser/ Skróty/ Abreviações/ Abrevieri/ Сокращения/ Skratky/ Okrajšave/ Förkortningar/ Kisaltmalar

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NC

NC means non-compliant and relates to radial balloon expansion which remains minimal between nominal and rated burst pressure

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Description

The Pantera LEO Fast-Exchange PTCA catheter (hereinafter Pantera LEO PTCA catheter or the device) is designed for dilatation of stenotic segments in coronary arteries and for post-dilatation of coronary stents. The dilatation balloon is designed to inflate to a known diameter and length at recommended inflation pressures (refer to Compliance Chart on page 2). There are two radiopaque markers, one at the proximal end and one at the distal end of the cylindrical part of the balloon, to facilitate fluoroscopic visualization and positioning of the catheter. The proximal end of the catheter has a single female Luer-Lock port to inflate/deflate the balloon. The working length of Pantera LEO PTCA catheter is 145 cm. The proximal shaft of the catheter is a hypotube and the distal portion of the catheter is a flexible shaft. The surface of the catheter is coated as described in Table 1 below:

Segment (mm)	Coating type
Balloon [Ø 2.0 – 5.0]	Hydrophobic
Distal shaft	Hydrophilic
Proximal shaft	Hydrophobic

The guide wire lumen starts at the catheter tip and ends at the guide wire exit port, 34 cm from the distal end. The Pantera LEO PTCA catheter is compatible with guide wires of 0.014" [0.36 mm] diameter and guiding catheters with an inner diameter \geq 0.056" [1.42 mm]. To indicate when the catheter tip exits from the guiding catheter, shaft exit markers are located on the hypotube 92 cm (brachial technique) and 102 cm (femoral technique) from the distal end of the Pantera LEO PTCA catheter. The hypotube fastener on the hub is designed to help manage and secure the proximal shaft/hypotube section of the PTCA catheter when it is gathered for storage or disposal. Do not attempt to hold the distal shaft using the hypotube fastener.

How Supplied

Sterile. Non-pyrogenic. Device is sterilized with ethylene oxide.

Contents

- One Pantera LEO PTCA catheter and one Compliance Data Card in a sealed, peel open pouch.
- One Instructions for Use Manual.

Storage

Store protected from sunlight and keep dry. Store between 15°C [59°F] and 25°C [77°F]. Short-term storage temperature excursions between 10°C [50°F] and 40°C [104°F] are allowed up to 48 hours.

Indications

The Pantera LEO PTCA catheter is indicated for balloon dilatation of the stenotic portion of a coronary artery or bypass graft stenosis for the purpose of improving myocardial perfusion and for post-dilatation of coronary stents.

Contraindications

All general contraindications for percutaneous transluminal coronary angioplasty (PTCA) are contraindications for this device. Contraindications for this device and PTCA in general are:

- Lesions that cannot be reached or treated with the system.
 - Target vessel unprotected left main coronary artery.
 - Coronary artery spasm in the absence of a significant stenosis.
- Furthermore, all procedure-related contraindications as described in the national and international guidelines of the respective medical associations apply.

Warnings

- This device is designed and intended for single use only. DO NOT resterilize and/or reuse. Reuse of single-use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness, or death of the patient. Cleaning, disinfection and sterilization may compromise essential material and design characteristics leading to device failure. BIOTRONIK will not be responsible for any direct, incidental or consequential damages resulting from resterilisation or reuse.

- Before use, visually inspect the packaging. If the packaging is damaged or opened, or if any information is obscured, DO NOT use the device.
- DO NOT use the device after the Use by date indicated on the label.
- DO NOT expose the catheter to organic solvents, e.g., alcohol.
- Use only an appropriate balloon inflation medium (e.g., 50:50 mixture by volume of contrast medium and saline). Never use air or any gaseous medium to inflate the balloon.
- When the catheter is in the body, it should be manipulated under high-quality fluoroscopy.
- DO NOT advance or retract the catheter unless the balloon is fully deflated under vacuum. If strong resistance is met during manipulation, stop the procedure and determine the cause of the resistance before proceeding. Advancement by force may result in damage to the vessel and/or laceration or separation of the guide wire or dilatation catheter. This may necessitate recovery of fragments.
- DO NOT exceed the rated burst pressure (RBP) indicated in the Compliance Chart. Use of a pressure-monitoring inflation device is mandatory to prevent over-pressurization.
- To reduce the potential for vessel damage, the inflated diameter of the balloon should approximate the original diameter of the vessel proximal and distal to the stenosis.

Precautions

- Prior to the procedure, the device should be visually examined to verify its functionality and to ensure that its size is suitable for the specific procedure for which it is to be used.
- This device should only be used in medical facilities by physicians thoroughly trained and experienced in performing vascular interventions, (including cases of life-threatening complications).
- Exercise care during handling to reduce the possibility of accidental breakage, bending or kinking of the catheter shaft.
- Ensure that the guide wire exit port remains inside the guiding catheter at all times. Guide wire exit port is indicated on the label.
- Appropriate anticoagulant, antiplatelet and vasodilatation therapy should be administered to the patient.
- Use only guide wires with a 0.014" [0.36 mm] diameter. When using two guide wires, care should be taken when introducing, torquing and removing one or both guide wires to avoid entanglement or spiralling. Caution should be applied not to damage the Pantera LEO PTCA catheter with the second guide wire during the advancement and the retrieval.
- Use guiding catheters with an inner diameter \geq 0.056" [1.42 mm].
- The hypotube fastener is designed to secure the proximal shaft/hypotube section of the PTCA catheter. The distal shaft should not be held by the hypotube fastener.
- DO NOT use the Pantera LEO PTCA catheter if a vacuum cannot be held, as this indicates a system failure.

Potential Adverse Events/Complications

Possible adverse events associated with PTCA procedure include, but are not limited to:

- Cardiac events: Myocardial infarction or ischemia, abrupt closure of coronary target vessel, restenosis of treated vessel, cardiogenic shock, angina, tamponade, perforation or dissection of coronary artery or aorta, cardiac perforation, pericardial effusion, aneurysm formation.
- Arrhythmia events: Ventricular tachycardia, ventricular fibrillation, atrial fibrillation, bradycardia.
- Balloon catheter events: Failure to reach or cross the lesion, inflation difficulties, rupture or pinhole of the balloon, deflation difficulties, withdrawal difficulties, embolization of catheter material.
- Respiratory events: Acute pulmonary edema, congestive heart failure, respiratory insufficiency or failure.
- Vascular events: Access site hematoma, hypotension/hypertension, pseudoaneurysm, arteriovenous fistula formation, retroperitoneal hematoma, vessel dissection, rupture or perforation, intimal tear, vessel recoil and restenosis, thrombosis or occlusion, vasospasm, peripheral ischemia, distal embolization (air, tissue debris, thrombus).
- Neurologic events: Permanent (stroke) or reversible (transient ischemic attack) neurologic event, femoral nerve injury, peripheral nerve injury.
- Bleeding events: Access site bleeding or hemorrhage, hemorrhage requiring transfusion or other treatment.
- Allergic reactions to contrast media, antiplatelets, anticoagulants.
- Death.
- Renal disorder.
- Infections.

Directions for Use

Select an appropriate Pantera LEO PTCA catheter for the target vessel. The inflation diameter of the balloon must not exceed the diameter of the coronary artery or bypass graft proximal or distal to the stenosis. If the stenosis cannot be crossed with the desired dilatation catheter, a smaller diameter PTCA catheter can be used to predilate the lesion.

PTCA Catheter Preparation

1. Take the protection ring with the catheter out of the package and place it onto a sterile field.
2. Gently pull out the catheter from the protection ring.
3. Remove the balloon protector and the stylet.

Flushing of the Guide Wire Lumen

4. Connect a syringe containing sterile saline to an appropriately sized flushing needle.
5. Apply the needle into the distal tip of the catheter and flush the guide wire lumen.
6. Remove the syringe and the flushing needle.

Remove Air from the Catheter

7. Connect a 3-way stopcock to the Luer-Lock of the catheter.
8. Prepare and remove air from the inflation device according to the manufacturer's recommendations and instructions.
9. Attach the inflation device containing 3 ml of balloon inflation medium to the stopcock.

Warning: Use only an appropriate balloon inflation medium (e.g., 50:50 mixture by volume of contrast medium and saline). Never use air or any gaseous medium to inflate the balloon.

10. Open the stopcock so that the fluid path between the catheter and the inflation device is established.
11. Create a vacuum on the inflation device and remove air from the catheter for at least 30 seconds.
12. **Caution:** DO NOT use the Pantera LEO PTCA catheter if a vacuum cannot be held, as this indicates a system failure.
13. Close the stopcock so that the fluid path to the catheter is closed and evacuate all air from the inflation device through the stopcock.
14. Repeat steps 10-12 if necessary, to ensure that any air contained in the balloon and inflation lumen is removed. Release the inflation/deflation barrel to normal pressure.
15. Open the stopcock and set the system aside for use.

Insertion Technique

15. Attach a hemostatic valve to the Luer port of the guiding catheter positioned within the vasculature.
16. Position the guide wire under fluoroscopy, in accordance with PTCA techniques.
17. Thread the proximal end of the guide wire, into the distal tip of the PTCA catheter until it exits at the guide wire exit port 34 cm from the distal tip.
18. Carefully insert the PTCA catheter through the hemostatic valve.
19. Advance the PTCA catheter through the guiding catheter using fluoroscopic guidance to determine when the catheter tip approaches the distal tip of the guiding catheter.

Note: The two proximal exit markers may be used to approximate when the PTCA catheter has reached the distal end of the guiding catheter.

20. Advance the PTCA catheter into the target vessel and following the guide wire toward the lesion.

Note: The radiopaque balloon markers facilitate balloon positioning within the lesion.

Balloon Inflation

21. Inflate the balloon to dilate the lesion using standard PTCA techniques.

Warning: DO NOT exceed RBP.

Warning: To reduce the potential for vessel damage, the inflated diameter of the balloon must not exceed the original diameter of the vessel proximal and distal to the lesion.

22. If a significant stenosis persists, inflate the balloon again increasing pressure gradually until the lesion fails to improve.

23. After each inflation, assess the distal coronary blood flow by arteriography through the guiding catheter.

Balloon Deflation

24. Deflate the balloon in accordance with the standard PTCA procedures. Apply negative pressure to the balloon for at least 40 seconds before carefully pulling the PTCA catheter out of the target vessel.

25. Pull the completely deflated PTCA catheter into the guiding catheter.

Catheter Removal/Exchange Procedure

24. Loosen the hemostatic valve.

27. Hold the guide wire and hemostatic valve in one hand, while grasping the PTCA catheter shaft in the other hand.

28. Maintain the guide wire position in the coronary artery by keeping the guide wire stationary, and begin pulling the dilatation catheter out of the guiding catheter.

Note: Monitoring the guide wire position under angiographic guidance is highly recommended during the exchange.

29. Pull on the catheter until the guide wire exit port is reached. Carefully remove the flexible, distal portion of the PTCA catheter off the guide wire while maintaining the guide wire position across the lesion. Close the hemostatic valve.

Note: In case of difficulties during catheter removal, remove the entire system at once, i.e., the guiding catheter, the guide wire and the PTCA catheter simultaneously.

30. Completely remove the PTCA catheter from the guide wire.

31. If required, prepare and introduce the next PTCA catheter to be used, as previously described or following the indications of the manufacturer.

Note: Inspect the PTCA catheter integrity immediately upon removal from the patient.

32. After use, the device and packaging may be contaminated with infectious substances (e.g. blood). Therefore, dispose of the device and packaging in accordance with hospital, administrative and/or local government policy.

Warranty/Liability

This product and each of its components (hereinafter product) were designed, manufactured, tested and packaged with all reasonable care. However, since BIOTRONIK does not have control over the conditions under which the product is used, the contents of this Instructions For Use (IFU) are to be considered as an integral part of this disclaimer for cases when a disturbance of the product's intended function may occur for various reasons.

BIOTRONIK does not guarantee that the following events will not occur:

- Product malfunctions or failures
- Patient's immune response to the product
- Medical complications during the use of the product or as a consequence of the product being in contact with the patient's body
- BIOTRONIK does not assume any liability for:
 - The use of the product that is not in accordance with the stated intended use/indication, contraindications, warnings, precautions and the directions for use of this IFU
 - Modification to the original product
 - Events which could not have been foreseen at the time of product delivery using the available levels of science and technology
 - Events originating from other BIOTRONIK products or products not from BIOTRONIK and
 - Force majeure events including, without being limited to, natural disasters.

The above provisions shall be without prejudice to any disclaimer and/or limitation of liability agreed separately with the customer to the extent permitted by applicable laws.

Intellectual Property

The following list includes trademarks or registered trademarks of the BIOTRONIK Group of Companies in the EU, the United States, China, Japan, and possibly in other countries: Pantera, BIOTRONIK. All other trademarks are the property of their respective owners.

Safety and Clinical Performance

The Summary of Safety and Clinical Performance (SSCP) can be found at <https://ec.europa.eu/tools/eudamed> using the Basic UDI-DI: 76401304BUDI0000267.

Clinical Benefits

Post-dilatation of insufficiently or poorly apposed stents and dilation of stenosis when needed results in improved clinical outcomes.

Reporting of Serious Incidents

Report any serious incident that has occurred with the device to the manufacturer and to the competent authority of the country in which your hospital is established. In case of a device malfunction, stop using the device immediately and return it to the manufacturer.

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All rights reserved. Specifications subject to
BIOTRONIK AG
© 2021 BIOTRONIK AG