# Pantera Pro

<b>English</b>	Coronary Dilatation Catheter	Instructions for Use	6
<b>Deutsch</b>	Koronardilatationskatheter	Gebrauchsanweisung	7
<b>Français</b>	Cathéter de dilatation coronaire	Mode d'emploi	9
<b>Italiano</b>	Catetere per dilatazione coronarica	Istruzioni per l'uso	10
<b>Español</b>	Catéter de dilatación coronaria	Instrucciones de uso	12
<b>Български</b>	Коронарен дилатационен катетър	Инструкции за употреба	13
<b>Hrvatski</b>	Koronarni dilatacijski kateter	Upute za uporabu	15
<b>Česky</b>	Koronární dilatační katetr	Pokyny k použití	16
<b>Dansk</b>	Koronart dilatationskateter	Brugsanvisning	17
<b>Nederlands</b>	Dilatatiekatheter voor coronairarteriën	Gebruiksaanwijzing	19
<b>Suomi</b>	Sepelvaltimon laajennuskatetri	Käyttöohjeet	20
<b>Ελληνικά</b>	Καθετήρας διαστολής στεφανιαίων	Οδηγίες χρήσης	22
<b>Magyar</b>	Szív koszorúér-tágító katéter	Használati utasítás	24
<b>Latviešu</b>	Koronārais dilatācijas katetrs	Lietošanas norādījumi	25
<b>Lietuvių</b>	Vainikinių arterijų plečiamasis kateteris	Naudojimo instrukcija	26
<b>Norsk</b>	Koronardilatasjonskateter	Bruksanvisning	28
<b>Polski</b>	Cewnik rozszerzający do naczyń wieńcowych	Instrukcja użycia	29
<b>Português</b>	Cateter de dilatação coronária	Instruções de utilização	31
<b>Română</b>	Cateter de dilatare coronariană	Instrucțiuni de utilizare	32
<b>Русский</b>	Катетер для баллонного расширения коронарных артерий	Инструкция по применению	34
<b>Slovenčina</b>	Koronární dilatačný katéter	Návod na použitie	35
<b>Slovenščina</b>	Koronarni dilatacijski kateter	Navodila za uporabo	37
<b>Svenska</b>	Koronar dilatationskateter	Bruksanvisning	38
<b>Türkçe</b>	Koroner Dilatasyon Kateteri	Kullanma Talimatı	40

**Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.**

STERILE EO

**BIO BIOTRONIK**  
excellence for life

Compliance Chart/Elastizitätstabelle/Tableau de compliance/Tabella dei valori di cedevolezza del palloncino/Tabla de distensibilidad/Таблица за съвместимост/Grafikon rastezljivosti/Tabulka kompatibility/Komplianskort/Compliantietabel/venyvystautilukko/Διάγραμμα ενδοτικότητας/Tabúlekonysági táblázat/Tabula atbilstības noteikšana/Atitikties lentelė/Komplianstabel/ Tabela podatności/ Tabela de elasticidade/Tabel de conformitate/Таблица растяжимости/Tabulka podajnosti/Preglednica skladnosti/elasticitetstabel/ Esneklik Tablosu									
Inflation Pressure (IP)		Balloon Diameter (mm) [B Ø]							
	atm	[kPa]	1.25	1.5	2.0	2.5	3.0	3.5	4.0
	6	[608]	1.22	1.46	1.97	2.41	3.00	3.53	3.83
<b>NP</b>	7	[709]	1.24	1.49	2.01	2.49	3.08	3.62	3.95
	8	[811]	1.26	1.53	2.05	2.58	3.17	3.72	4.07
	9	[912]	1.27	1.56	2.09	2.66	3.24	3.80	4.18
	10	[1013]	1.30	1.60	2.12	2.72	3.31	3.87	4.28
	11	[1115]	1.31	1.63	2.16	2.79	3.37	3.93	4.36
	12	[1216]	1.33	1.66	2.18	2.84	3.42	3.97	4.43
	13	[1317]	1.35	1.69	2.21	2.88	3.45	4.02	4.49
<b>RBP</b>	14	[1419]	1.37	1.72	2.23	2.93	3.50	4.06	4.55
	15	[1520]	1.39	1.74	2.26	2.97	3.54	4.10	4.60
<b>NP</b>	In vitro testing has shown that all balloons reached their nominal size at given Nominal Pressure.								
<b>RBP</b>	In vitro testing has shown that with 95% confidence, 99.9% of the balloons will not burst at or below the given Rated Burst Pressure. DO NOT exceed RBP.								

## Abbreviations/Abkürzungen/Abréviations/Abbreviazioni/Abreviaturas/Абревиатури/Kratice/Zkratky/Forkortelser/Afkortingen/Lyhenteet/Συμπίναξ/Rövidítések/Saisinjumi/Santrumpos/Forkortelser/Skróty/Abreviaturas/Abrevieri/Сокращения/Skratky/Okrajšave/Förkortningar/Kisaltmalar

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Pantera Pro	Balloon Length (mm) [BL]					
Balloon Ø (mm) [B Ø]	6	10	15	20	25	30
1.25	X	X	X	X		
1.5	X	X	X	X		
2.0		X	X	X	X	X
2.5		X	X	X	X	X
3.0		X	X	X	X	X
3.5		X	X	X	X	X
4.0		X	X	X	X	X
Unavailable Product Sizes [UPS]						
Hydrophilic Coating of Balloon [HCB 1]						
Hydrophobic Coating of Balloon [HCB 2]						

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## Description

The Pantera Pro Coronary Dilatation Catheter [hereinafter Pantera Pro catheter] is intended to be used for dilatation of stenotic segments in coronary arteries or bypass grafts. The dilatation balloon is designed to inflate to a known diameter and length at recommended inflation pressures [refer to Compliance Chart]. There are two radiopaque markers, one at the proximal end and one at the distal end of the cylindrical part of the balloon, to facilitate fluoroscopic visualization and positioning of the catheter. The 1.25 mm and 1.5 mm balloon diameter versions have a single marker positioned in the middle of the cylindrical part of the balloon. The proximal end of the catheter has a single female Luer port to inflate/deflate the balloon. The working length of the Pantera Pro is 140 cm. The proximal shaft of the catheter is a hypotube and the distal portion is a flexible shaft. The surface of the catheter is coated as described in Table below:

Coating surfaces of the Pantera Pro catheter sizes:	
Segment (mm)	Coating type
Balloon [Ø 1.25 – 2.00]	Hydrophilic
Balloon [Ø 2.50 – 4.0]	Hydrophobic
Distal shaft	Hydrophilic
Proximal shaft	Hydrophobic

The guide wire lumen starts at the catheter tip and ends at the guide wire exit port, 29 cm from the distal end. The Pantera Pro is compatible with guide wires of 0.014" [0.36 mm] diameter and guiding catheters with an inner diameter of  $\geq 0.070$ " [1.78 mm; 6F]. In addition, two Pantera Pro catheters with balloon diameters not larger than 3.5 mm are compatible for a kissing balloon procedure using a guiding catheter with an inner diameter of  $\geq 0.070$ " [1.78 mm; 6F]. To indicate when the catheter tip exits from the guiding catheter, shaft exit markers are located on the hypotube 92 cm [brachial technique] and 102 cm [femoral technique] from the distal end of the Pantera Pro.

## How Supplied

Sterile. Non-pyrogenic. Device is sterilized with ethylene oxide. DO NOT use if the package is opened or damaged, or if any information provided is obscured or damaged.

## Contents

- One (1) Pantera Pro in a sealed, peel-open pouch.
- One (1) Instructions for Use Manual.
- One (1) Compliance Chart.

## Storage

- Store in a dark, dry, cool place
- Keep dry
- Keep away from sunlight

## Indications

Pantera Pro is indicated for balloon dilatation of the stenotic portion of a coronary artery or bypass graft stenosis for the purpose of improving myocardial perfusion. The Pantera Pro [balloon diameter 2.0 – 4.0 mm] is also indicated for post-delivery expansion of balloon expandable stents.

**Note:** Post-deployment stent expansion testing was performed on the bench with the PRO-Kinetic Energy coronary stent system. All stents should be deployed in accordance with the manufacturer's indications and instructions for use.

## Contraindications

All general contraindications for percutaneous transluminal coronary angioplasty (PTCA) are contraindications for this device.

Contraindications for this device and for PTCA in general are:

- Lesions that cannot be reached or treated with the system.
- Target vessel unprotected left main coronary artery.
- Coronary artery spasm in the absence of a significant stenosis.

Furthermore, all procedure-related contraindications as described in the national and international guidelines of the respective medical associations apply.

## Warnings

- This device is designed and intended for single use only. DO NOT sterilize and/or reuse. Reuse of single-use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of the patient. Cleaning,

disinfection and sterilization may compromise essential material and design characteristics leading to device failure. BIOTRONIK will not be responsible for any direct, incidental or consequential damages resulting from resterilization or reuse.

- DO NOT use if the package is opened or damaged or if any information provided is obscured or damaged.
- DO NOT expose the catheter to organic solvents, e.g., alcohol.
- To reduce potential for air embolus into the vessel, use only an appropriate balloon inflation medium [e.g., 50:50 mixture by volume of contrast medium and saline]. Never use air or any gaseous medium to inflate the balloon.
- DO NOT use the catheter if a vacuum cannot be held, as this indicates a system failure.
- When the catheter is in the body, it should be manipulated only under high-quality fluoroscopy.
- DO NOT torque the catheter more than one (1) full turn.
- DO NOT advance or retract the catheter unless the balloon is fully deflated under vacuum. If resistance is met during manipulation, stop the procedure and determine the cause of the resistance before proceeding. Advancement by force may result in damage to the vessel and/or laceration or separation of the guide wire or dilatation catheter. This may necessitate recovery of fragments.
- DO NOT exceed the Rated Burst Pressure (RBP). The rated burst pressure is based on the results of in vitro testing. The use of a pressure monitoring device is mandatory to prevent overpressurization.
- To reduce the potential for vessel damage, the inflated diameter of the balloon must not exceed the original diameter of the vessel proximal and distal to the lesion. When a kissing balloon procedure is performed make sure that the combined sizes of two inflated balloons do not exceed the diameter of the vessel.
- PTCA in patients who are not acceptable candidates for coronary artery bypass graft surgery warrants careful consideration, including possible hemodynamic support during PTCA, as treatment of this patient population carries special risk.

## Precautions

- Prior to procedure, the catheter should be examined to verify functionality and ensure that its size and shape is suitable for the specific procedure for which it is to be used. Special attention in selecting suitable catheter sizes should be applied when kissing balloon procedure is performed.
- Only physicians thoroughly trained and educated in the performance of percutaneous transluminal coronary angioplasty (PTCA) should use the catheter system.
- Exercise care during handling to reduce the possibility of accidental breakage, bending or kinking of the catheter shaft.
- Use prior to the "use by" date specified on the packaging.
- Ensure that the guide wire exit port remains inside the guiding catheter at all times. Guide wire exit port is indicated on the label.
- Appropriate anticoagulant, antiplatelet and vasodilatation therapy should be administered to the patient.
- Use only guide wires with a 0.014" [0.36 mm] diameter. When using two guide wires, care should be taken when introducing, torquing and removing one or both guide wires to avoid entanglement or spiralling. Caution should be applied not to damage the Pantera Pro with the second guide wire during the advancement and the retrieval.
- Use guiding catheters with an inner diameter  $\geq 0.056$ " [1.42 mm; 5F]. For kissing balloon procedures, use two Pantera Pro catheters with balloon diameters not larger than 3.5 mm within a guiding catheter with an inner diameter  $\geq 0.070$ " [1.78 mm; 6F]. Advance or retrieve the catheters sequentially [NOT simultaneously]. Care should be taken when Pantera Pro is used in combination with another type of dilatation balloon catheter.
- If the surface of the distal portion of the Pantera Pro becomes dry, the coating should be reactivated by immersion in sterile saline.
- After use, disposal of product and product packaging must be executed in accordance with hospital, administrative and/or local government policy.

## Adverse Effects

Possible adverse events associated with the use of Pantera Pro when used as intended include, but are not limited to, the following:

- Cardiac events: Acute myocardial infarction or ischemia, acute closure of coronary target vessel, restenosis of treated vessel, cardiogenic shock, unstable angina, tamponade, perforation or dissection of coronary artery or aorta, cardiac perforation, emergency cardiac surgery, pericardial effusion, aneurysm formation, palpitations.
- Arrhythmia events: Ventricular tachycardia, ventricular fibrillation, atrial fibrillation, bradycardia.
- Balloon catheter events: Failure to reach or cross the lesion, inflation difficulties, rupture or pinhole of the balloon, deflation difficulties, withdrawal difficulties, embolization of catheter material.

- Respiratory events: Acute pulmonary edema, congestive heart failure, respiratory insufficiency or failure.
- Vascular events: Access site hematoma, hypotension/ hypertension, pseudoaneurysm, arteriovenous fistula formation, retroperitoneal hematoma, vessel dissection, rupture or perforation, intimal tear, vessel recoil and re-stenosis, thrombosis or occlusion, vasospasm, peripheral ischemia, distal embolization [air, tissue debris, thrombus].
- Neurological events: Permanent (stroke) or reversible (TIA) neurologic event, femoral nerve injury, peripheral nerve injury.
- Bleeding events: Access site bleeding or hemorrhage, hemorrhage requiring transfusion or other treatment.
- Allergic reactions to contrast media, antiplatelets, anticoagulants.
- Infections.
- Death.

## Directions for Use

Select an appropriate Pantera Pro for the target vessel.

**Warning:** The inflated diameter of the balloon must not exceed the original diameter of the vessel proximal and distal to the lesion.

If the stenosis cannot be crossed with the desired dilatation catheter, a smaller diameter dilatation catheter can be used to pre-dilate the lesion.

**Caution:** For kissing balloon procedures two Pantera Pro catheters with balloon diameters not larger than 3.5 mm can be advanced sequentially (NOT simultaneously) within a guiding catheter with an inner diameter  $\geq 0.070$ " [1.78 mm; 6F].

## Dilatation catheter preparation

1. Take the protection ring with the catheter out of the package and place onto a sterile field.
2. Gently pull out the catheter from the protection ring.
3. Carefully remove the balloon protector together with the transportation wire.

**Warning:** DO NOT use the catheter if unusual resistance is felt during the removal of the transportation wire and/or balloon protector.

## Pre-flush guide wire lumen

4. Connect a syringe containing sterile saline to an appropriately sized flushing needle.
5. Apply the needle into the distal tip of the catheter and flush the guide wire lumen.
6. Remove the syringe and the flushing needle.

## Purge air from the catheter

7. Connect a 3-way stopcock to the Luer port of the catheter.
8. Prepare and remove air from a 20 ml capacity inflation/deflation device according to manufacturer's recommendations and instructions.
9. Attach the inflation/deflation device containing 3 ml of balloon inflation medium to the stopcock.

**Warning:** Use only an appropriate balloon inflation medium [e. g., 50:50 mixture by volume of contrast medium and saline].

**Never use air or any gaseous medium to inflate the balloon.**

10. Open the stopcock so that the fluid path between the catheter and the inflation/deflation device is established.
11. Pull the plunger of the inflation/deflation device and aspirate air from the catheter for at least 30 seconds.

**Warning:** DO NOT use the catheter if a vacuum cannot be held, as this indicates a system failure.

12. Close the stopcock so that the fluid path to the catheter is closed and evacuate all air from the inflation/deflation device through the stopcock.
13. Repeat steps 10–12 if necessary, to ensure air contained in the balloon and inflation lumen are removed. Release the inflation/deflation barrel to normal pressure.
14. Open the stopcock and set the system aside for use.

## Insertion technique

15. Attach a hemostatic valve to the Luer port of the guiding catheter positioned within the vasculature.
16. Position the guide wire under fluoroscopy, in accordance with PTCA techniques.
17. Immerse the distal portion of the catheter in sterile saline to activate the hydrophilic coating. Ensure the guiding catheter is also thoroughly flushed with sterile saline.
18. Back load the proximal end of the guide wire into the distal tip of the dilatation catheter until it exits at the guide wire exit port 29 cm from the distal tip.
19. Carefully insert the dilatation catheter through the hemostatic valve.
20. Advance the dilatation catheter through the guiding catheter using fluoroscopic guidance to determine when the catheter tip approaches the distal tip of the guiding catheter.

**Note:** The two proximal exit markers may be used to approximate when the dilatation catheter has reached the distal end of the guiding catheter.

21. Advance the dilatation catheter into the target vessel and follow the guide wire towards the lesion.
- Note:** The radiopaque marker(s) facilitate(s) balloon positioning within the lesion.
- Warning:** If resistance is met during manipulation, stop the procedure and determine the cause of the resistance before proceeding. Advancement by force may result in damage to the vessel and/or laceration or separation of the guide wire or the dilatation catheter. This may necessitate recovery of fragments.

## Balloon inflation

22. Inflate the balloon to dilate the lesion using standard PTCA techniques.
- Warning:** Do not exceed a total of 10 inflations.
- Warning:** DO NOT exceed Rated Burst Pressure (RBP).
- Warning:** To reduce the potential for vessel damage, the inflated diameter of the balloon must not exceed the original diameter of the vessel proximal and distal to the lesion.
23. If a significant stenosis persists, reinflate the balloon with increasing pressure until the stenosis improves.
24. After each inflation, assess the distal coronary blood flow by arteriography through the guiding catheter.

## Balloon deflation

25. Deflate the balloon in accordance with the standard PTCA procedures. Apply negative pressure to the balloon for at least 30 seconds before carefully pulling the dilatation catheter out of the target vessel.
26. Pull the completely deflated dilatation catheter into the guiding catheter.

## Catheter removal/Exchange procedure

27. Loosen the hemostatic valve.
28. Hold the guide wire and hemostatic valve in one hand, while grasping the dilatation catheter shaft in the other hand.
29. Maintain the guide wire position in the coronary artery by keeping the guide wire stationary, and begin pulling the dilatation catheter out of the guiding catheter while monitoring the guide wire position under fluoroscopy.
30. Pull on the catheter until the guide wire exit port is reached. Carefully remove the flexible, distal portion of the dilatation catheter off the guide wire while maintaining the guide wire position across the lesion. Close the hemostatic valve.
- Caution:** In case of difficulties during catheter removal, remove the entire system at once, i.e., the guiding catheter, the guide wire and the dilatation catheter simultaneously.
31. Completely remove the dilatation catheter from the guide wire.
32. If required, prepare and introduce the next dilatation catheter to be used, as previously described or following the indications of the manufacturer.
- Note:** Inspect the dilatation catheter integrity immediately upon removal from the patient.
- Caution:** Disposal of packaging materials can be handled per normal industrial/local standards. The used device shall be handled according to hospital procedure.

## Warranty/Liability

This product and each of its components (hereinafter product) were designed, manufactured, tested and packaged with all reasonable care. However, since BIOTRONIK does not have control over the conditions under which the product is used, the contents of this Instructions For Use (IFU) are to be considered as an integral part of this disclaimer for cases when a disturbance of the product's intended function may occur for various reasons.

BIOTRONIK does not guarantee that the following events will not occur:

- Product malfunctions or failures
- Patient's immune response to the product
- Medical complications during the use of the product or as a consequence of the product being in contact with the patient's body.

BIOTRONIK shall assume no liability for:

- The use of the product that is not in accordance with the stated intended use/indication, contraindications, warnings, precautions and the directions for use of this IFU
- Modification to the original product
- Causes which could not have been foreseen at the time of product delivery using the available levels of science and technology
- Causes originating from other BIOTRONIK products or products not from BIOTRONIK and force major events which include, but not limited to, natural disasters.

The above provisions shall be without prejudice to any disclaimer and/or limitation of liability agreed separately with the customer to the extent permitted by applicable laws.

For US only:

BIOTRONIK disclaims all warranties, expressed or implied regarding the product, including but not limited to, any warranty of merchantability

or fitness for a particular purpose of the product. BIOTRONIK is not liable for any direct, incidental or consequential damages or expenses caused by any use of the product whether the claim is based on any theory of liability, including negligence.

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All rights reserved. Specifications subject to
BIOTRONIK AG
© BIOTRONIK AG