


Vascular Intervention // **Coronary**
Rapid Exchange Microcatheter

Micro Rx

Transforming guidewire support and
lesion crossing with Rx convenience



Expand the reach of your workhorse guidewire¹



Simple rapid exchange support to cross
tortuous vessels & highly stenosed lesions*¹



Ultra-low profile** to facilitate lesion entry
and penetration¹

Micro Rx


Rapid Exchange Microcatheter

Indicated to support a guidewire during access of coronary and/or peripheral vasculature.^a

Vascular
Intervention
Coronary



Transforming guidewire support and lesion crossing with Rx convenience



Technical Data

Rapid exchange microcatheter

Tip design	Tapered soft tip
Lesion entry profile	0.44 mm / 0.017" / 1.3F
Shaft design	Reinforced shaft
Tip exit markers	95 cm and 105 cm
Proximal exit markers	110 cm and 120 cm
Coating	Hydrophilic, distal 12 cm
Guidewire compatibility	0.014" / 0.36 mm

Ordering Information

Ordering number	IMDS article number	Guiding catheter compatibility	Usable length
483054	RX135	≥ 5F	135 cm

^aIndication as per IFU.

1. IMDS data on file. *Non CTO **compared to FineCross

Micro Rx is a trademark of IMDS. Distributed by BIOTRONIK in selected countries.

MICRO Rx

Manufacturer:
IMDS Operations B.V.
Ceintuurbaan Noord 150
9301 NZ Roden,
The Netherlands
Tel +31 (0)50 8200230
info@imds.nl

Distributor:
BIOTRONIK, Inc.
6024 Jean Road
Lake Oswego, OR 97035-5369
Tel: (800) 547-0394 (24-hour)
Fax: (800) 291-0470
www.biotronik.com



©2023 BIOTRONIK, Inc. All rights reserved.
SS254 December 2023

